



**Connect**

- 1** Connect your USB cable between the COM port on the Vapor Sorption Analyzer (VSA) and the USB port on your computer
- 2** Select the proper COM port in the toolbar
- 3** Click the Connect button in the toolbar

## Setup an Isotherm Test

- 1** Once your VSA is connected, click the Wizard button in the toolbar
- 2** Follow the Wizard's instructions and fill in the parameters of your test
- 3** After completing the Wizard, follow the instructions displayed on the VSA to start test

**Download and Export**

- 1** Once your VSA is connected and a test has been run, select Download in the toolbar
- 2** From the File menu select "Export..." to create an Excel compatible table
- 3** To view your data in table format, click the Table button in the toolbar
- 4** To view your data in chart format, click the Chart button in the toolbar

# **Access Analysis Tools**



- 1** From the Tools menu dropdown at the top of the screen, select the desired tool.
- 2** Many tools require the use of isotherm information and use a .isotherm file. This file is created using the Generate Isotherm Tool, which appears first in the list of tools. These files will then be utilized in other tools.
- 3** Enter the required information for the selected tool and then activate to obtain results.

---

While the Moisture Analysis Toolkit does not have a hardcopy manual, extensive help files exist to guide you through the various features and processes.

---

---

**Contact Us:** 509-332-5601  
support@aqualab.com

**Check for Updates:** [aqualab.com/softwareupdates](http://aqualab.com/softwareupdates)

---

#####

11-12-14



Copyright ©2014 Decagon Devices, Inc.  
2365 NE Hopkins Court, Pullman, WA 99163

Printed in USA