




METER

AQUALAB® TDL/TDL 2

BENCHTOP WATER ACTIVITY METER

**Installation Qualification and Operational
Qualification (IQ/OQ) Protocols and Instructions**





METER

TABLE OF CONTENTS

1. Introduction.....	3
2. Installation Qualification (IQ).....	4
2.1 Initial Qualification and Requalification	4
2.2 Equipment Identification.....	4
2.3 Receiving and Unpacking.....	5
2.4 Visual Inspection.....	6
2.5 Environmental Conditions.....	7
2.6 Power-Up Test.....	8
3. Operational Qualification (OQ)	9
3.1 Hardware Testing.....	9
3.2 System Information.....	9
3.2.1 Sensor Verification.....	9
3.2.2 Equipment Familiarization.....	11
Appendix A: Training Record.....	12
Appendix B: Deficiencies and Corrective Actions.....	13



METER

1. Introduction

This qualification protocol is solely intended to be used with new or relocated AQUALAB® TDL or TDL 2 Water Activity Meters. It is written to assist the user in validation of predetermined specifications. The use of this document does not replace the need for the AQUALAB TDL User Manual.

Information within the user manual is required to complete this IQ/OQ Protocol. If the manual has been misplaced, copies can be obtained from the manufacturer or downloaded from metergroup.com/TDL-support. Qualification of instrumentation is a formal process of documenting that an instrument is fit for its intended use and that it is kept maintained and calibrated.

Responsibilities

The instrument qualification carried out onsite is the sole responsibility of the instrument owner/user. However, METER Group, Inc. USA supports customers in performing the qualification by providing the instrument qualification, dedicated documentation, and offering a qualification service. In this regard, the following responsibilities are defined:

Performance of Qualification

Execution of the instrument qualification and entire qualification of the installed system covered in this document is performed by authorized METER Group, Inc. USA service specialists when ordered from a customer.

Review and Final Qualification Approval

Final approval for the qualification must be completed after review of the qualification documentation filled out during performance of the qualification procedures (IQ/OQ protocols). The customer representative then signs the approved form.

Installation Qualification (IQ)

Installation qualification is documented proof that the instrument was received as designed and specified by the manufacturer, that it is properly installed in the selected environment, and that this environment is suitable for the operation and use of the instrument. The IQ section, therefore, describes and documents the instrument installation in the predetermined environment. Further, the IQ verifies and ensures that all ordered parts and documentation are in place and that all supplied items are in working order and condition.

Operational Qualification (OQ)

The operational qualification serves as proof that the equipment operates as designed and intended, as well as fulfills acceptance criteria defined and stated in the Operational Qualification documentation. These criteria are defined and are based on the equipment technical specifications of the manufacturer.

Performance Qualification (PQ)

Performance qualification is documented proof that an instrument consistently performs according to the specifications appropriate for its routine use. Monitoring of equipment during routine operation is essential for ensuring that the ongoing performance is within specifications. The performance qualification, execution, and frequency are the sole responsibility of the user. Performance validation should be designed to meet the specifications and accuracy for a given application.

Equipment Familiarization and Operator Training Records

All equipment users are to be instructed in basic operation, functionality, instrument parameters, as well as on basic hardware features of the installed system, including routine maintenance and cleaning procedures. Please contact METER Group, Inc. USA to learn about available training and seminars. Authorized support specialists perform the qualification services offered by METER Group, Inc. USA.

Please identify which instrument this IQ/OQ applies to by marking either AQUALAB TDL or TDL 2.

AQUALAB TDL

AQUALAB TDL 2

The IQ/OQ will not commence without acknowledgement of METER Group, Inc. USA *Terms and Conditions*, signature, and date below. Please fill out this form as completely as possible. Please refer to metergroup.com/terms-conditions for details.

I agree to the METER Group, Inc. USA Service terms, Payment terms, and Terms and Conditions.

Signature _____

Date _____

Name/Title _____



METER

2. Installation Qualification (IQ)

The Installation Qualification (IQ) protocols described in this section are dedicated to initial qualification and/or to requalification. IQ tests should be performed in the following instances:

1. When the system is installed
2. When the system is moved to a new location
3. Prior to running Operational Qualification (OQ) tests

This section describes the procedure for receiving, unpacking, and installing an AQUALAB TDL/TDL 2 Water Activity Meter.

2.1 Initial Qualification and Requalification

Fill out the company information below after unpacking the AQUALAB TDL/TDL 2 Water Activity Meter and corresponding accessories.

Company Information

Company _____

Department _____

Street Address _____

City _____

State _____

Zip _____

Contact Information

First Name _____

Last Name _____

Email _____

Phone _____

2.2 Equipment Identification

Fill out this section after unpacking the AQUALAB TDL/TDL 2 Water Activity Meter and corresponding accessories.

Manufacturer METER Group, Inc. USA _____

Model _____

Serial Number _____

METER Group, Inc. USA
Authorized Representative

Signature _____

Date _____

Name/Title _____

Initials _____



METER

2.3 Receiving and Unpacking

Fill out the information below, verifying the condition of the external packaging (make note if any damage during shipment in a way might have damaged the internal package content).

External packaging conditions Satisfactory Not Satisfactory

Notes and remarks about receiving and packaging.

Checking Contents

Compare shipment list with supplied items to ensure completeness of order.

AQUALAB TDL Benchtop Unit	<input type="checkbox"/> Complete	<input type="checkbox"/> Not Complete	<input type="checkbox"/> NA
Quick Start Guide	<input type="checkbox"/> Complete	<input type="checkbox"/> Not complete	<input type="checkbox"/> NA
User Manual	<input type="checkbox"/> Complete	<input type="checkbox"/> Not Complete	<input type="checkbox"/> NA
Certificate of Calibration	<input type="checkbox"/> Complete	<input type="checkbox"/> Not Complete	<input type="checkbox"/> NA
Trial Verification Standards	<input type="checkbox"/> Complete	<input type="checkbox"/> Not Complete	<input type="checkbox"/> NA
SDS Documents	<input type="checkbox"/> Complete	<input type="checkbox"/> Not Complete	<input type="checkbox"/> NA
Power Cable	<input type="checkbox"/> Complete	<input type="checkbox"/> Not Complete	<input type="checkbox"/> NA
USB Cable	<input type="checkbox"/> Complete	<input type="checkbox"/> Not Complete	<input type="checkbox"/> NA
Cleaning Kit	<input type="checkbox"/> Complete	<input type="checkbox"/> Not Complete	<input type="checkbox"/> NA
AQUALAB Qualification Kit	<input type="checkbox"/> Complete	<input type="checkbox"/> Not Complete	<input type="checkbox"/> NA
LDPE Sample Cups and Lids	<input type="checkbox"/> Complete	<input type="checkbox"/> Not Complete	<input type="checkbox"/> NA

Notes and remarks from checking contents.



METER

Order Completeness

All parts were received as ordered and the delivery is complete. Yes No

Any parts that were missing at the time of this supply verification and reported as Not Complete must be delivered to complete the shipment. Any parts marked as Not Complete must be indicated below and reviewed with the customer. Minor parts or accessories that do not impact the installation or qualification procedure or the functionality of the meter can be accepted, if agreed upon by the customer, to complete the remainder of the IQ/OQ process. Completed can be marked once the item has been received.

		Item	Date	Initials
<input type="checkbox"/> Accepted	<input type="checkbox"/> Completed			
<input type="checkbox"/> Accepted	<input type="checkbox"/> Completed			
<input type="checkbox"/> Accepted	<input type="checkbox"/> Completed			
<input type="checkbox"/> Accepted	<input type="checkbox"/> Completed			
<input type="checkbox"/> Accepted	<input type="checkbox"/> Completed			
<input type="checkbox"/> Accepted	<input type="checkbox"/> Completed			

Notes and remarks about order completeness.

2.4 Visual Inspection

After unpacking, verify that there is no physical damage to the meter, cables, and accessories. Note all observed damage in the notes and remarks section below. Minor defects that do not affect functionality can be marked as accepted, if approved by the customer.

AQUALAB TDL / TDL 2 Meter	<input type="checkbox"/> Satisfactory	<input type="checkbox"/> Not Satisfactory	<input type="checkbox"/> Accepted
Documentation	<input type="checkbox"/> Satisfactory	<input type="checkbox"/> Not Satisfactory	<input type="checkbox"/> Accepted
Other	<input type="checkbox"/> Satisfactory	<input type="checkbox"/> Not Satisfactory	<input type="checkbox"/> Accepted

Notes and remarks about the visual inspection.

Qualified by:

Signature	Date
Name/Title	Initials



METER

2.5 Environmental Conditions

The AQUALAB TDL / TDL 2 Water Activity Meter should be installed on a level surface in a location where the temperature remains stable. This location should be well away from air-conditioner and heater vents, open windows, outside doors, or other items that may cause rapid temperature fluctuation or vibration.

Location	<input type="checkbox"/> Satisfactory	<input type="checkbox"/> Not Satisfactory	<input type="checkbox"/> Accepted
Adequate Power	<input type="checkbox"/> Satisfactory	<input type="checkbox"/> Not Satisfactory	<input type="checkbox"/> Accepted
Stable Surface	<input type="checkbox"/> Satisfactory	<input type="checkbox"/> Not Satisfactory	<input type="checkbox"/> Accepted
Temperature	<input type="checkbox"/> Satisfactory	<input type="checkbox"/> Not Satisfactory	<input type="checkbox"/> Accepted

Notes and remarks about the visual inspection.

Qualified by:

Signature _____

Date _____

Name/Title _____

Initials _____



METER

2.6 Power-Up Test

After finding a suitable location for the AQUALAB TDL/TDL 2, plug the power cord into the back of the unit and the other end into a standard AC outlet.

AQUALAB TDL /TDL 2 is powered upon switching on. Yes No

Notes and remarks about the power-up test.

Completeness of Installation Qualification (IQ)

IQ was completed and documented according to METER guidelines. Yes No

Notes and remarks about the completeness of the IQ.

METER Group, Inc. USA Authorized Representative

Signature	Date
<hr/>	<hr/>
Name/Title	Initials
<hr/>	<hr/>

IQ was reviewed by the representative of the system owner.

Qualified by:

Signature	Date
<hr/>	<hr/>
Name/Title	Initials
<hr/>	<hr/>



METER

3. Operational Qualification (OQ)

This section contains descriptions of tests that are to be executed for the Operational Qualification (OQ) of the AQUALAB TDL/TDL 2 Water Activity Meter to prove proper operation.

3.1 Hardware Testing

The display is functional Yes No
 The keypad is functional Yes No

Notes and remarks about the visual inspection.

Qualified by:

Signature	Date
_____	_____
Name/Title	Initials
_____	_____

3.2 System Information

Manufacturer METER Group, Inc. USA	Firmware Version
_____	_____
Model	Serial Number
_____	_____

3.2.1 Sensor Verification

The AQUALAB TDL/TDL 2 Water Activity Meter utilizes a tunable diode laser to determine water activity. The performance of the laser is verified by measuring specially prepared calibration standards that have a specific molality and water activity. Performance Verification Standards in four water activity levels are used for qualification: 0.250 a_w , 0.500 a_w , 0.760 a_w , 1.000 a_w . The AQUALAB TDL will read each standard within $\pm 0.005 a_w$ of the stated value. To measure the water activity of the standards, follow the instructions for taking a reading in the user manual.

Standard @ 25 °C	Lot	$a_w \pm 0.005$	Tunable Diode Laser °C
13.41m LiCl	0.250	_____	_____
	0.250	_____	_____
8.57m LiCl	0.500	_____	_____
	0.500	_____	_____
6.0m NaCl	0.760	_____	_____
	0.760	_____	_____
Steam Distilled H ₂ O	1.000	_____	_____
	1.000	_____	_____

Water Activity Standards are within specification. Yes No



METER

Notes and remarks about the verification standards measurement inspection.

If verification standards are out of specification, clean the AQUALAB TDL/TDL 2 and follow the procedure in the user manual to perform a linear offset. Repeat the verification process with fresh standards.

Standard at °C	Lot	$a_w \pm 0.005$	Tunable Diode Laser °C
13.41m LiCl	0.250	<hr/>	<hr/>
	0.250	<hr/>	<hr/>
8.57m LiCl	0.500	<hr/>	<hr/>
	0.500	<hr/>	<hr/>
6.0m NaCl	0.760	<hr/>	<hr/>
	0.760	<hr/>	<hr/>
Steam Distilled H ₂ O	1.000	<hr/>	<hr/>
	1.000	<hr/>	<hr/>

Water Activity Standards are within specification. Yes No

Notes and remarks about the verification standards measurement inspection.

Qualified by:

Signature	Date
Name/Title	Initials



METER

3.2.2 Equipment Familiarization

This section ensures that the instrument operators receive appropriate equipment training to ensure proper operation, maintenance, and generation of results with the AQUALAB TDL/TDL 2 Water Activity Meter. Product familiarization covers instruction on basic operations, functionality, and meter features, and routine maintenance, including cleaning procedures.

Equipment familiarization and training was completed for the AQUALAB TDL. Yes No

Qualified by:

Signature	Date
_____	_____
Name/Title	Initials
_____	_____

Completeness of Operational Qualification (OQ)

OQ was completed and documented according to METER guidelines. IQ Requalification

Qualification met the Vendor Acceptance Criteria Yes No

NOTE: If any deficiencies are found, fill out the instructions for a corrective action on [page 13](#) of this document.

METER Group, Inc. USA Authorized Representative

Signature	Date
_____	_____
Name/Title	Initials
_____	_____

Operational Qualification was reviewed by the representative of the system owner.

Qualified by:

Signature	Date
_____	_____
Name/Title	Initials
_____	_____



METER

Appendix A: Training Record

This training record is for instruction in basic operation, functionality, meter parameters, as well as on basic hardware features of the installed system, including routine maintenance and cleaning procedures. Please contact METER Group, Inc. USA to learn about available training and seminars.

Authorized support specialists perform the qualification services offered by METER Group, Inc. USA.

Name	Date
_____	_____
Name	Date
_____	_____
Name	Date
_____	_____
Name	Date
_____	_____
Name	Date
_____	_____
Name	Date
_____	_____
Name	Date
_____	_____
Name	Date
_____	_____
Name	Date
_____	_____
Name	Date
_____	_____
Name	Date
_____	_____

**METER Group, Inc. USA
Authorized Representative**

Signature	Date
_____	_____
Name/Title	Initials
_____	_____

Operational Qualification was reviewed by the representative of the system owner.

Qualified by:

Signature	Date
_____	_____
Name/Title	Initials
_____	_____



METER

Appendix B: Deficiencies and Corrective Actions

If any deficiencies were found, a corrective action will be opened. Once acceptable results are obtained, the results of the corrective action are accepted by checking the Accepted box below.

Deficiency 1

Explain the deficiency.

Corrective Action

The deficiency was reviewed, and corrective action was taken and accepted. Accepted

Deficiency 2

Explain the deficiency.

Corrective Action

The deficiency was reviewed, and corrective action was taken and accepted. Accepted

METER Group, Inc. USA
Authorized Representative

Signature _____

Date _____

Name/Title _____

Initials _____