

Document Title: <b>Slick, Pawkit Chinese</b>		Part # and Rev. <b>18103</b>	
		Release Date:	
Rev.	Description	Revision By	Date
Date -- Time	Added description from other slicks	Allison	6/10/16

**Production Filename:**

<http://publications.decagon.com/Marketing/Description%20Files/18103 Slick Pawkit CN.pdf>

**Printing Process:** Digital

**Size:** 11 inches wide, 8.5 inches tall

**Ink:** 4/4 Color Both Sides

**Paper:** 100 lbs. book gloss coated

**Page Count:** 2

**Proof type:** soft proof

**Single hole drilling:** none

**AquaLab Pawkit 便携式水分活度仪**

**AQUA LAB**  
BY DECAGON

**最小的水分活度仪**

- 低维护, 传感器和读数一体化
- 快速检测, 测量时间小于5分钟
- 易读数, 水分活度和样品温度共同显示
- 电池可用3年, 无需外接电源

- 超轻便, 易于携带
- 一年质保
- 校准简便
- 包括: 手提箱, 样品杯, 标准溶液, 快速启动指南以及手册



**专为卫生监督员、食品科学家以及食品、药品和化妆品等质检人员的快速检查而设计**

**长10 cm重115 g的水分活度仪**

AquaLab Pawkit是一款真正的便携式水分活度仪, 长度只有10cm, 重量只有115g, 和手机一样轻巧。Pawkit在5分钟内测量水分活度, 并达到±0.02 a<sub>w</sub>的精度。内置锂电池可使用三年, 最适于质检人员随身携带, 随时随地测量。没有专业实验室的生产企业也可用它来做快速检测。



实验室外测量的理想选择