

Document Title: Description, Serial, Series 4, AquaLab Label		Part # 13491	
		Release Date: 07/06/07	
Rev.	Description	Revision By	Date
01	Updated Logo, Font	DDH	3/31/11
02	Updated Logo, Font, barcode	DDH	9/14/12
03	Add UL	NJR	5/20/16

Production Filename:

http://publications.decagon.com/Vendors/13491_Series4Serial%20Label_Print.pdf

Dimensions: 2.3125 inch wide, 1.5 inch tall

Material: DPM Poly White Sheet 2.0 mil

Colors: Black ink only

Printer: Eclipse Screen Printing

Special Notes: Illustration is not to scale ** reference only **

Serial number: Start serial # at

Barcode: Add Barcode and SN# under Logo - please add 2 separate smaller barcodes

